

Crock meter Digital

Purpose:

To find out 'colour fastness to rubbing' of any kind of fabric /yarn/ garments.

Can used as per both ISO and AATCC standards

Features:

- Automatic rotations are made with help of electric motor to avoid human errors.
- All componets are made with mild steel and chrome plated for long life and trouble- free operations.
- Test specimen is made to rub against standard adjacent fabric under standard test conditions.



Model No.	: CT/CMD 104	Ringer	: 2nos
Counter	: Fully Automatic Digital Counter	Length of Rubbing:	104 as per standard
Applicable Force	: 9 newton		
Finger Dia			
Specimen Holder	: 1 no Heavy Metal		