

Hot Air Oven Cum Incubator

Purpose:

Used to perform Perspiration and Oven dry test.

Features:

Used in perspiration & oven dry tests.

Digital temperature controller incorporated.

18" X 18" X 18" inside chamber dimension.

Outside mild steel – beautifully painted and fully stainless-steel inner chamber for rust free operations.

Supplied with adjustable stainless-steel trays to hold samples.

Circulating fan is provided for uniformity.

Accuracy maintained by using energy regulator.

Load indicator indicates heater status.

1 blower to maintain uniform temperature inside chamber.

Related Standards

AATCC 15, 106, 107, ISO 105 E02, E01, E04

