

## ICI Pilling Tester 2 Station

### **Purpose:**

Used to replicate pilling in fabrics under standard test conditions as per ICI method.

### **Features:**

Model : CT/PT2

Pilling tester is used to replicate pilling in fabrics under standard test conditions as per ICI method.

Pilling tester with digital LCD controls. A unit comes with 2 pilling boxes with standard cork lining fitted inside and 5-digit pre-set counter.

The boxes are rotated about their central axis as per standard conditions.

Inside, size of each box is (235 X 235 X 235) mm without cork lining thickness of 3.2 mm.

Pilling photographs – ASTM or IWS are available optionally.

Option of 4 station also available. With 1 No. Specimen Mounting Jig and 1 No. of Specimen Marking Template. Per box 4 Rubber Test Tube.

Related Standards

ISO 5470, ASTM D4966

