

Lab Shaker

Purpose:

Used in laboratories for continuous shaking of solutions / liquid

Features:

Model No : CTRS

With speed regulation system to change the level of shaking.

Up 12,16,24,32 beakers can be kept together for shaking. multiple samples can be tested at a time.

Specially designed stainless-steel beaker holders give flexibility rust free operations to use glass beakers without damage.

Digital timer with memory Back and display (Pre-Set Counter).

Digital RPM meter for controlling and indicating RPM of the top rotating unit.

30±1 RPM, 60±1 RPM, 90±1 RPM can be achieved.

Smooth precision-engineered components for exceptional performance.

Used for fibre composition tests.

Power Supply 230 VAC, 50Hz

