

# Pilling Assessment Viewer

## **Purpose:**

To assess and grade pilled samples against standards

## **Features:**

Offers standard conditions to view effect of pilling.

Suitable for all standards where the assessment of pilling on fabrics is necessary, whether against control fabrics or photographs.

It uses D65 Artificial Daylight (with CWF optionally available).

Suitable for: Martindale Pilling, ICI Pilling, ICI Snagging, Random Tumble Pilling, Brush/Sponge Pilling.

## **Related Standards**

ISO 12945, ASTM D3512, ASTM D4970, ASTM D5362, IWSTM 152, BS 5811

