

Print Durability Tester – Durawash

Purpose:

Used for repeating washing/ colour fastness tests or specific garment or fabric print durability tests, pleat retention, waddings and quilted waddings

Features:

The Durawash simulates the conditions necessary to meet leading chain stores standard garment and fabric durability tests.

The process control systems provide excellent accuracy and repeatability.

Multi range timer with end of cycle audible alarm automatic agitator shut off.

Display resolution 0.1-degree C.

System control and calibration accuracy better than +/-1 °C.

Process tank construction in high grade stainless-steel.

Industrial quality motor specifications.

Model No.	: CT/ PD /560
Spin Type	: As per Standard
Rotational Spee	: As per Standard
Heating Elements	: 4000W
TemperatureRange	: 0-299°C Pre Selectable Type (PT100)
Timer Range	: 0-99.00 Mins (Minutes/Seconds Pre Selectable Type)

